Tiinuse test piimast

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IDEXX Milk Pregnancy Test

Best practices for the collection and handling of milk samples

The quality of samples taken for any diagnostic procedure is extremely important. Correct identification and a truly representative sample are critical for obtaining optimal diagnostic test results. Below is a protocol for the manual collection of milk samples for chemical-based pregnancy testing. Use bovine milk samples from ≥ 28 days postbreeding and ≥ 60 days postcalving. For optimal timing of sample collection in other species, refer to the test's user instructions.

- Materials for sampling**
- 1. Clean vials or tubes 30 mL–50 mL capacity (Standard herd testing tubes will suffice).
 - 2. Preservative (e.g., tablet containing bronopol) pre-filled in the vial or tube
 - 3. Wet wipes or towel
 - 4. Permanent marker or labels for identifying sample vials/tubes

IDEXX Milk Pregnancy Test
Best Practices for the Collection of Milk Samples

The quality of samples taken for any diagnostic procedure is extremely important. Correct identification and a truly representative sample are critical for obtaining optimal diagnostic test results. Below is a protocol for the manual collection of milk samples for chemical-based pregnancy testing.

Materials for Sampling

- 1. Sterile vials or tubes 30–50mL capacity (standard herd testing tubes will suffice).
- 2. Preservative e.g. Bronopol or micro-tablets (if available)
- 3. Wet-wipes or towel
- 4. Permanent marker or labels for identifying sample vials/tubes

Sampling Technique

- 1. Label tube prior to sampling. Include date and individual cow ID. It is essential to ensure that the correct tube is used for the correct animal. Failure to do so may result in an incorrect result being assigned to the cow.
- 2. Using the wet-wipes or towel, remove dirt and detritus from the teats. Ensure udder is clean and dry prior to starting sample collection.
- 3. Samples may be taken from individual quarters or as a composite sample of all four quarters.
- 4. Discard a few streams of milk and remove the cap from the tube or vial. Angle tube at 45 degrees to the teat end.
- 5. Three-quarters fill the tube or vial from one or multiple quarters and replace the cap.
- 6. Place samples in a refrigerator prior to shipping to the laboratory.
- 7. Samples should be well packaged and kept cool during shipping (especially if preservative has not been used).

Sample Quality and Handling at the Laboratory

- 1. Bronopol or a similar preservative may be used to maintain sample quality—this will not affect test results.
- 2. The assay can be run on whole or skim milk samples and samples that have been exposed to heat treatment (during milk component analysis).
- 3. Poor quality samples may compromise the accuracy of test results. Prior to testing, milk samples will be checked to ensure that they are not soured or separated and are free from contamination.

Contact IDEXX Technical Services if you have any questions. Within the United States and Canada, call 1-800-548-9977. Outside of the United States and Canada, call 1-327-356-4899.

Sampling technique

- 1. Label tube prior to sampling. Include date and individual cow ID. It is essential to ensure that the correct tube is used for the correct animal. Failure to do so may result in an incorrect status being assigned to the cow.
- 2. Using the wet wipes or towel, remove dirt and detritus from the teats. Ensure udder is clean and dry prior to starting sample collection.
- 3. Samples may be taken from individual quarters or as a composite sample of all four quarters.
- 4. Discard a few streams of milk and remove the cap from the tube or vial. Angle tube at 45° to the teat end.
- 5. Fill the tube or vial three-quarters full from one or multiple quarters and replace the cap.
- 6. Place samples in a refrigerator prior to shipping to the laboratory.
- 7. Samples should be well-packaged and kept cool during shipping (especially if preservative has not been used).

Sample quality and handling at the laboratory

- Bronopol or a similar preservative may be used to maintain sample quality—this will not affect test results.
- The assay can be run on whole or skim milk samples and samples that have been exposed to heat treatment (during milk component analysis, usually up to 20 minutes at 40°C/105°F).
- Poor-quality samples may compromise the accuracy of test results. Prior to testing, milk samples will be checked to ensure that they have not soured or separated and are free from contamination.
- Do not heat milk samples in microwave as this procedure will destroy the pregnancy-associated glycoproteins (PAGs).

Sample type	Sample storage conditions and time frames		
	Room temperature	Refrigerator (2°C–8°C)	Freezer (-20°C)
Native milk samples	Less than a few hours	3 days	Indefinite
Preserved milk samples (bronopol or other preservative)	3–5 days	1 month	Indefinite

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Test with Confidence.™

Lihtne ja usaldusväärne

TRUSTED BY FARMERS AROUND THE GLOBE

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




Tiinuse testide kasv (proovid)

Tiinuse test piimast Taanis 2018 (2013-2018)

RYK provides milk recording for dairy cows. RYK are collecting annually about 5.5 million milk samples, serves 3.000 dairy farmers, and have a turnover of 16 million euro. We have 65 employees, and offices in Aarhus, Sorø, Holstebro, and Vojens.

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1. Choose PCR or Pregnancy test

2. Choose methode

3. Choose criteria for pointet out cows methode B

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

Graph shows SEENOREST growth from January 2017 to June 2018

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Küsimused??

Test with Confidence.™

IDEXX

Viited

1. Paul Fricke, Steve Stewart, Paul Rapnicki, Steve Eicker, and Michael Overton Pregnant vs. Open: Getting Cows Pregnant and the Money it Makes
2. De Vries et al, Economics of improved reproductive performance in dairy cattle (Publication AN 156), Gainesville FL, University of Florida Institute of Food and Agriculture Science; 2005
3. A. Ricci , P. D. Carvalho , M. C. Amundson , R. H. Fourdraine , L. Vincenti , and P. M. Fricke; Factors associated with pregnancy-associated glycoprotein (PAG) levels in plasma and milk of Holstein cows during early pregnancy and their effect on the accuracy of pregnancy diagnosis
4. Vasconcelos, J. L. M., R. W. Silcox, J. A. Lacerda, J. R. Pursley, and M. C. Wiltbank. 1997. Pregnancy rate, pregnancy loss, and response to heat stress after AI at 2 different times from ovulation in dairy cows. Biol. Reprod.